

ReSound Remote Control



Product description

The ReSound Remote Control is for clients who appreciate the ability to easily and discreetly control their hearing aids.

They can use it to change programs and adjust volume, without even looking. This is thanks to a simple keypad that features large buttons, which are concave, tactile and ergonomic. They have raised symbols in the centre, so they are easy to locate and glide a finger over. This is especially good for users who struggle with manual dexterity.

The ReSound Remote Control is small and handy and has a streamlined, elegant design. It fits in a handbag, a pocket or the palm of the hand like a keyless car key – and will easily blend in with other smart accessories on a desk or a dinner table.

TECHNICAL DATA

Dimensions	73 x 35 x 14 mm
Weight	22.5 g
Power supply	CR2032 Lithium coin cell battery
Wireless connections	2 hearing aids per Remote Control
Wireless frequency	2.4 GHz
Wireless range	Approximately arm's length distance
Communication speed	2 Mbit/s
Emission	Max. +4 dBm, in compliance with all relevant standards
ESD tolerance	+/- 4 kV contact. +/- 8 kV air.
Operating temperature	0 to +55 °C / +32 to +131 °F
Storage temperature	-20 to +60 °C / -4 to +140 °F

EASY TO CONNECT

- To pair with a hearing aid, press the pairing button on the back of the ReSound Remote Control and open and close the battery door of the hearing aid within 20 seconds.

EASY TO USE

- You can change programs in your hearing aid with a simple push on the program button.
- Your hearing aid volume can be adjusted with the + and – button.
- The standard CR2032 Lithium coin cell battery is widely available and easy to replace.

Worldwide Headquarters
GN Hearing A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: +45 4575 1111
resound.com

United Kingdom
GN Hearing UK Ltd.
Kirtlington Business Centre
Portway, Kirtlington
Oxon OX5 3JA
Tel.: +44 1869 352 800
resound.com

Australia
GN Hearing Australia Pty Ltd
Gate C, 19-25 Khartoum Road
Macquarie Technology Park
Macquarie Park NSW 2113
Tel.: (free) 1800 658 955
resound.com

New Zealand
GN Hearing New Zealand Limited
Ground Floor, North Entrance
4 Fred Thomas Drive
Takapuna, Auckland, 0622
Tel.: (free) 0800 900 126
resound.com

Singapore
GN Hearing Pte. Ltd.
2 Kallang Avenue
#07-19 CT HUB
Singapore - 339407
Tel.: +65 6320 9388
resound.com

CVR no. 55082715